



Technical Data

3M™ HB Quat Disinfectant - Ready-to-Use

Description

3M™ HB Quat Disinfectant is ready to dispense using any cloth, mop, sponge or sprayer.

This is a broad spectrum hospital disinfectant for use on hard, non-porous surfaces such as floors, walls, desks or countertops, tables and chairs, telephones, bed frames and cabinets to clean, disinfect and deodorize. They are effective in 300 ppm hard water (calculated as CaCO₃) in the presence of organic soil (5% blood serum). The ready-to-use solution contains 708 ppm quat.

Applications

Use on washable, hard, non-porous environmental surfaces in:

- Hospitals, surgi-centers, nursing homes, doctor and dentist offices, including patient care rooms, newborn nurseries
- Schools, health clubs, spas, restrooms, and shower rooms
- Veterinary clinics and dairies
- Restaurants, bars and food processing plants

Packaging

3M™ HB Quat Disinfectant Ready-to-Use Rinse Free

- 5-gallon bottle
- 50-gallon Drum
- 250-gallon Tote
- 4 50-gallon Drums per pallet

Directions for Use

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

Pour 3M™ HB Quat Disinfectant into a properly labeled sprayer bottle or mop bucket for application as desired.

BACTERICIDAL STABILITY OF USE-DILUTION:

Tests confirm that this ready to use product in the presence of 5% soil load, remains effective against *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Salmonella enterica* for up to one year when stored in a sealed container at room temperature.

If the ready use product becomes visibly dirty or contaminated, the use product must be discarded and replaced with fresh product. Always use clean, properly labeled dry containers when using the product. Bactericidal stability of the product does not apply to open containers such as buckets or pails. Product in open containers must be prepared daily or more often if the solution becomes visibly dirty, diluted, or contaminated.

NOTE:

Refer to product Safety Data Sheets for specific physical properties, health hazard, first aid and precautionary information.

- 18-1459-9 (RTU)

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that: (1) introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes, but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body.

This product may be used to pre-clean or decontaminate critical or semi-critical devices prior to sterilization or high-level disinfection.

Disinfection/Fungicidal*Virucidal*

Directions:

Apply to hard inanimate, non-porous surfaces thoroughly wetting surfaces with a cloth, mop, sponge or sprayer. For heavily soiled areas, a preliminary cleaning is required. For sprayer applications use a coarse spray device. Spray 6-8 inches from surface and rub with brush, sponge or cloth. Do not breathe spray. Use to disinfect hard, non-porous surfaces. Treated surfaces must remain wet for 10 minutes. Product used in open container must be prepared at least daily or earlier when it becomes diluted or soiled.

General Deodorization:

To deodorize, pour into a properly labeled bottle or mop bucket. Excess material must be wiped up or allowed to air dry.

When used as directed, this product will deodorize surfaces in restroom and toilet areas, behind and under sinks and counters, garbage cans and garbage storage areas, and other places where bacterial growth can cause bad odors.

KILLS HIV-1, HCV and HBV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD OR BODY FLUIDS

in health care settings or other settings in which there is an expected likelihood of soiling with blood or body fluids, and in which the surfaces or objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (AIDS Virus), Hepatitis C Virus (HCV), and Hepatitis B Virus (HBV).

Special instructions for cleaning and decontamination of HIV-1, HCV & HBV on surfaces/objects soiled with blood or body fluids:

Personal Protection:

Specific barrier protection items to be used when handling items soiled with blood or body fluids are disposable protective gloves, gowns, masks and eye protection.

Cleaning Procedure:

Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of this product.

Disposal of infectious material:

Blood and other body fluids, cleaning materials and clothing must be autoclaved or disposed of according to Federal, State and local regulations for infectious waste disposal.

Contact Time:

Leave surfaces wet for 30 seconds for HIV-1 and 10 minutes for HCV and HBV. The contact time for the viruses, fungi and bacteria listed on this label is 10 minutes except for Polio virus Type 1 (chat strain) which is 30 minutes.

Disinfection and Virucidal*

Performance:

HB Quat Disinfectant Cleaner 25 is bactericidal, virucidal*, and fungicidal on hard inanimate surfaces modified in the presence of 5% organic serum with a 10 minute contact time unless otherwise noted below and found to be effective against the following microorganisms on hard, non-porous environmental surfaces:

Isolates From AIDS Patients

- *Aspergillus niger*
- *Candida albicans*
- *Cryptococcus neoformans*
- *Pseudomonas aeruginosa*
- *Staphylococcus aureus*
- *Streptococcus pneumoniae*

Gram Positive Clinical Isolates

- *Enterococcus faecalis*
- *Micrococcus luteus*
- *Staphylococcus aureus*
- *Staphylococcus epidermidis*
- *Staphylococcus saprophyticus*
- *Streptococcus haemolyticus*
- *Streptococcus pyogenes*

Gram Negative Clinical Isolates

- *Acinetobacter calcoaceticus* var. *anitratus*
- *Acinetobacter calcoaceticus* var. *lwoffii*
- *Bordetella bronchiseptica*
- *Burkholderia cepacia*
- *Brevundimonas dimunita*
- *Enterobacter agglomerans*
- *Enterobacter cloacae*
- *Enterobacter gergoviae*
- *Enterobacter liquefaciens*
- *Escherichia coli* (Urinary)
- *Escherichia coli* (Wound)
- *Flavobacterium meningosepticum*
- *Hafnia alvei*
- *Klebsiella oxytoca*
- *Klebsiella pneumoniae*
- *Morganella morganii*
- *Proteus mirabilis*
- *Proteus vulgaris*
- *Pseudomonas aeruginosa*
- *Pseudomonas fluorescens*
- *Pseudomonas pseudomallei*
- *Pseudomonas putida*
- *Pseudomonas stutzeri*
- *Serratia marcescens*
- *Sphingomonas paucimobilis*

Other Bacteria

- *Actinobacillus pleuropneumoniae*
- *Actinomyces pyogenes*
- *Bacillus cereus*
- *Bacteroides fragilis*
- *Corynebacterium ammoniagenes* (*Brevibacterium ammoniagenes*)
- *Bordetella bronchiseptica*
- *Burkholderia pickettii*
- *Campylobacter jejuni*
- *Chryseomonas luteola*
- *Corynebacterium pseudotuberculosis*
- *Enterobacter aerogenes*
- *Enterococcus faecalis*
- *Enterococcus faecium*
- *Enterococcus hirae*
- *Escherichia coli*
- *Escherichia coli* strain O157:H7
- *Escherichia vulneris*
- *Haemophilus influenzae*
- *Klebsiella pneumoniae*
- *Listeria monocytogenes*
- *Pasteurella haemolytica*
- *Pseudomonas aeruginosa*
- *Rhodococcus equi*
- *Salmonella enterica*
- *Salmonella schottmuelleri*
- *Salmonella typhi*
- *Shigella dysenteriae*



- Staphylococcus aureus
- Staphylococcus auricularis
- Staphylococcus capitis
- Staphylococcus hominis
- Staphylococcus simulans
- Stenotrophomonas maltophilia
- Streptococcus equi var. equi
- Streptococcus equi var. zooepidermicus
- Streptococcus pneumoniae (PRSP)
- Streptococcus pyogenes
- Streptococcus salivarius
- Yersinia enterocolitica

Pathogenic Fungi

- Trichophyton mentagrophytes

Environmental Fungi

- Aspergillus candidus
- Aspergillus niger
- Penicillium chermesinum
- Penicillium oxalicum
- Penicillium spinulosum
- Ulocladium sp.

Antibiotic Resistant Gram Negative Bacteria

- Pseudomonas aeruginosa (Sulfa, Cefatoxime, Nitrofurantoin, Tetracycline, Amikacin, Ampicillin, Cephalothin and Bactine Resistant)
- Escherichia coli (Ampicillin, Tetracycline, Penicillin and Sulfa Resistant)
- Klebsiella oxytoca (Ampicillin, Sulfanilimide and Tetracycline Resistant)
- Klebsiella pneumoniae type 1 (Ampicillin, Tetracycline, Cephalothin and Sulfa Resistant)
- Morganella morganii (Penicillin and Tetracycline Resistant)

- Enterobacter agglomerans (Ampicillin and Sulfanylimide Resistant)
- Salmonella choleraesuis (Antibiotic Resistant)
- Enterobacteriaceae with extended beta-lactamase resistance (Ampicillin and Piperacillin Resistant)

Antibiotic Resistant Gram Positive Bacteria

- Enterococcus faecalis (Vancomycin Resistant-VRE)
- Enterococcus faecium (Vancomycin Resistant-VRE)
- Staphylococcus aureus (Methicillin-MRSA, Community Associated Methicillin Resistant - CA-MRSA PVL Positive)
- Staphylococcus aureus (CA-MRSA Genotype USA 400)
- Staphylococcus aureus (Penicillin G, Penicillin, Ampicillin, Cefazolin, Cefatoxime, Chloramphenicol, Ciprofloxacin, Clindimycin, Erythromycin, Oxacillin, Rifampin, Tetracycline Resistant)
- Staphylococcus aureus (Vancomycin Resistant – VRSA)
- Staphylococcus aureus (Vancomycin Resistant Intermediate-VISA)
- Staphylococcus epidermidis (Ampicillin and Drug Resistant)

Human Viruses

- *Adenovirus type 2
 - *Cytomegalovirus
 - *HBV (Hepatitis B Virus)
 - *HBC (Hepatitis C Virus)
 - *Herpes Simplex Virus type 1
 - *Herpes Simplex Virus type 2
 - *HIV-1 (AIDS Virus)
 - *Human Coronavirus
 - *Influenza A/Brazil Virus
 - *Influenza A/Victoria (H3N2)
 - *Influenza A2-Asian Virus
 - *Influenza B Virus (Allen strain)
 - *Influenza C Virus (Taylor strain)
 - *Measles Virus
 - *Parainfluenza type 1
 - *Poliovirus type 1 (Chat strain)
- 30 minutes contact time

- *Respiratory Syncytial Virus
- *Rotavirus
- *Vaccinia Virus

Non-Human Viruses

- Avian Influenza/Turkey/Wisconsin Virus
- Canine Coronavirus
- Canine Distemper Virus
- Canine Herpesvirus
- Equine Herpesvirus
- Equine Influenza Virus
- Feline Calicivirus
- Norovirus
- Feline Infectious Peritonitis Virus
- Infectious Bovine Rhinotracheitis (IBR) Virus
- Newcastle Disease Virus
- Porcine Parvovirus
- Porcine Respiratory & Reproductive Syndrome Virus (PRRSV)
- Porcine Rotavirus
- Pseudorabies Virus
- Transmissible Gastroenteritis (TGE) Virus
- T1 bacteriophage
- T4 bacteriophage
- Vesicular Stomatitis Virus (VSV)
- Bovine Viral Diarrhea Virus (BVDV)
- Avian Influenza Virus (H5N1)
- Influenza A Virus (H1N1) (formerly called swine flu)

Product Specifications:

(Typical Values)

Active Ingredients:

Alkyl (60% C₁₄, 30% C₁₆, 5% C₁₂, 5% C₁₈) dimethyl benzyl ammonium chloride 13.238%

Alkyl (68% C₁₂, 32% C₁₄,) dimethyl ethylbenzyl ammonium chloride 13.238%

EPA Registration:

EPA Reg. No. 61178-5-10350

EPA Est. No. 10350-WI-001 (D), 75521-WI-001 (P) add 75521-WI-002 (P)



For order information contact your sales representative at 1-775-544-8662, or visit www.ds-gov.com/#products

Important Notice to User

Technical Information: The technical information, recommendations and other statements contained in this document are based upon tests or experience that 3M believes are reliable, but the accuracy or completeness of such information is not guaranteed.

Product Use: Many factors beyond 3M's control and uniquely within user's knowledge and control can affect the use and performance of a 3M product in a particular application. Given the variety of factors that can affect the use and performance of a 3M product, user is solely responsible for evaluating the 3M product and determining whether it is fit for a particular purpose and suitable for user's method of application.

Warranty, Limited Remedy and Disclaimer: 3M warrants that each 3M product meets the applicable product specification at the time the product is shipped and will be free from defects in material and manufacture for 90 days from the date of purchase from 3M's authorized distributor. **3M MAKES NO OTHER WARRANTIES OR CONDITIONS, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OR CONDITION OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR ANY IMPLIED WARRANTY OR CONDITION ARISING OUT OF A COURSE OF DEALING, CUSTOM, OR USAGE OF TRADE.** If a 3M product does not conform to this warranty, the sole and exclusive remedy is, at 3M's option, to replace the product or refund the purchase price.

Limitation of Liability: Except the extent prohibited by law, 3M will not be liable for any loss or damage arising from the 3M product, whether direct, indirect, special, incident or consequential, regardless of the legal or equitable theory asserted.



**Commercial
Solutions Division**
St. Paul, MN 55144-1000
1-800-852-9722
www.3M.com/facility